

EXPERIMENT 25: NEUTRALIZATION AND TITRATION

Equipment: 2 burets (marked acid and base), 2 small funnels (marked A and B), one 100mL or 150 mL beaker, glass stirring rod, iron ring, ring stand, buret holder, wire gauze, Bunsen burner
Materials: 20 mL 0.20M NaOH, phenolphthalein, 20mL HCl of unknown concentration, matches

In a titration a solution of known concentration is used to determine the concentration of an unknown solution.



- A.
 1. Fill a buret (labelled **BASE**; with a pinch clamp) with approximately 20 mL of a solution of NaOH of known molarity (approx. 0.20 M).
 2. Fill another buret (labelled **ACID**; with a glass stopcock) with 20 mL of a solution of HCl of unknown molarity. **Be sure the tips of the burets are filled.**
 3. Record the initial volumes in the table below. The meniscus (curved surface of the liquid) should always be read at eye level, and the readings should always be made at the bottom of the meniscus.

- B.
 4. Run 10mL of the acid into a 150 mL beaker.
 5. Add 1 drop of phenolphthalein indicator, and stir.
 6. Add the base drop-wise, and stir after each drop, until the solution remains pink.
 7. Read the meniscus on each buret and record the final volumes in the table below.

- C.
 8. Repeat the entire procedure in A and B at least once in order to test the precision of the method. Save the final solution from step B6 for part "D".

- D.
 9. Place the beaker containing a few mL of the neutralized solution on a wire gauze over a burner.
 10. Boil until all but a few drops of liquid have evaporated. Then shut off the burner, and allow the remainder to evaporate to dryness. **Caution:** Overheating may crack the beaker.
 11. What do you observe? Explain. _____

 12. Write an equation for the reaction between the acid and the base.

	TRIAL 1	TRIAL 2	TRIAL 3	TRIAL 4
Base: final reading				
Base: initial reading				
mL of base used				
Acid: final reading				
Acid: initial reading				
mL of acid used				

3. Calculate the molarity of the acid for each trial (to the correct number of significant figures).

$$M_{\text{base}} \times V_{\text{base used}} = M_{\text{acid}} \times V_{\text{acid used}}$$

Trial 1

Trial 2

Trial 3




Trial 4

Average acid Molarity: _____

4. Compare results obtained in the different titration trials.

5. List possible sources of experimental error.

6. Complete and balance the following equations:



7. How many milliliters of 0.50M HNO₃ would be needed to neutralize 15.0 mL of 1.00M KOH?

CONCLUSIONS

In a neutralization reaction an _____ and a _____ react to form a _____ and _____. Phenolphthalein, litmus, methyl orange, and methyl red are often used in titrations as acid-base _____.