

Lab Instructor _____
Date _____

LAB # _____ Plant and Animal Cells
Name _____
Period _____

Objective: To compare the structure of plant and animal cells

Use full sentences when answering all questions.

Pre-Lab

1. How big is a typical cell?
2. One millimeter contains 1000 microns. Convert the following. Show your work
 - a. How many millimeters are in 100 microns?
 - b. How many microns are in 5 millimeters?
 - c. Convert your height into microns. Show your work:

LAB

Materials


Chunks of red onion, compound light microscope, lens paper, slides, cover slips, prepared animal cell slides, forceps, transparent metric ruler, Lugol's solution, paper towels

Procedures and Observations

A. Observing characteristics in plant cells.


1. Take a piece of the inner layers of an onion chunk.
 - a. Use the forceps to carefully peel off a small piece of thin skin from the CONCAVE side of the layer. Make a temporary "wet mount" of this piece of the onion skin.
 - b. Examine the specimen using the 10X objective lens. Record all observations here.
2. Put a small piece of paper towel on one side of the cover slip. Carefully add some drops of the Lugol's solution on the opposite side. The paper towel should soak up extra liquid, hence pulling through the solution.
 - a. Observe the cells using the 10X objective lens. Draw the onion cells and label the cell wall, cytoplasm, vacuole and nucleus. Write the *total* magnification in the space provided.
 - b. How did the specimen change when Lugol's solution was added?
 - c. Carefully switch to the high power (40X) objective. Draw the onion cells and label the cell wall, cytoplasm, vacuole, cell membrane and nucleus. Write the *total* magnification in the space provided.

ONION EPIDERMIS (low power)



Magnification _____X

ONION EPIDERMIS (high power)



Magnification _____X

3. Remove the slide. Using the 10X objective, place the ruler across the widest part of the stage opening so that the millimeter edge faces you.
 - a. Focus the microscope so that the millimeter markings are clear.
 - b. Describe how the image of the ruler differs from its actual position.

Name _____

Period _____

4. Move the millimeter markings so that one of the lines is at the end of the field of diameter. Estimate the field diameter by noting the number of millimeter lines in view. What is the diameter of the field? Show all work.

In Millimeters: _____

In Micrometers: _____

5. Center the onion slide that you prepared over the stage opening.
 a. Use the 10X objective, count the number of cell lengths along the field diameter.
 b. Count the number of cell that fit width-wise across the field (from top to bottom).

Dimension	# of cells	Average length in mm	Average length in microns
Length (1 st count)			
Length (2 nd count)			
Width (1 st count)			
Width (2 nd count)			

6. Calculate an estimate of the average length and width of an onion cell in microns:

Length: _____


Width: _____

B. Observing animal cell characteristics.


1. Observe the prepared animal cell slide and draw and label the cytoplasm, vacuole cell membrane and nucleus.

ANIMAL CELL (low power)

ANIMAL CELL (high power)



Magnification _____X



Magnification _____X

2. Repeat Part A, step #5 using the prepared animal cell slide.

Dimension	# of cells	Average length in mm	Average length in microns
Length (1 st count)			
Length (2 nd count)			
Width (1 st count)			
Width (2 nd count)			

3. Calculate an estimate of the average length and width of an animal cell in microns:

Length: _____

Width: _____

4. Is it acceptable to use your estimations of cell size as a standard for use in scientific work? Why or why not?

5. How does your height in microns compare mathematically to the plant and animal cell lengths?

C. Compare and Contrast plant and animal cells.

Fill out the two charts using your observations in this lab and your prior knowledge.


ANIMAL CELL	PLANT CELL
HOW ARE THEY DIFFERENT?	
WITH REGARD TO	
ENERGY	
STRUCTURE & SIZE	
ORGANELLES	
MOBILITY	
FUNCTION	

Name _____ Period _____

ANIMAL CELL


PLANT CELL

HOW ARE THEY ALIKE?




A large empty rectangular box intended for the student to write their comparison of animal and plant cells.

PATTERNS OF SIGNIFICANT SIMILARITIES AND DIFFERENCES:



A large empty rectangular box intended for the student to write patterns of significant similarities and differences between animal and plant cells.

CONCLUSION OR INTERPRETATION:



A large empty rectangular box intended for the student to write their conclusion or interpretation of the comparison.