	NAME		PREF	SECTION SC	DATE
		EXPERIM	ENT 32: PR	RIMARY CELLS	
	Equipment: Materials:	150mL beakers, salt brid 5mL Cu ⁺² [aq], 5mL Pb ⁺² [aq],		Agar & KNO₃), sensitive vo u, Pb & Zn strips	oltmeter.
	In this expe	riment you will learn ho	w electricity ca	n be produced in a red	ox reaction.
A.	Place a strip of Cu in a beaker containing 1 M Cu (NO ₃) ₂ . In another beaker place a strip of Zn and 1 M Zn (NO ₃) ₂ solution. Connect the Cu and Zn strips to a voltmeter using single strand electric wire. 1. Why is there no voltage registered on the voltmeter?				
				Va	
В.	Make a salt bridge by soaking a strip of paper towel with a saturated KNO ₃ solution. Shake off the excess solution. Connect the two beakers using the salt bridge. Note which electrode is connected to the positive terminal of the voltmeter. 2. What voltage is registered on the voltmeter? 3. Using Reference Table N, determine which metal is oxidized, and which metal ion is				
	red	luced.			
	Write the half-reaction for each cell. 4. What is the theoretical voltage for a Cu/Cu (NO ₃) ₂ // Zn(NO ₃) ₂ /Zn cell?				
	E. Define and detically the toronto.				
	6. The 7. De	e electrode at which oxid fine reduction in terms of e electrode at which reduction	lation takes pla of electron trans	ice is called thesfer.	
C.		ts A and B using a Cu/			·
D.		nits, construct a Zn/Zn (cell. What voltage would		Diagram of E	lectro-chemical Cell
SU	MMARY QUE	STIONS			
	9. What ar	e practical applications arned in this experiment			
	10. List th	e possible sources of err ment.	or in this		
		e above cell, electrons flo electrode to the ode.			
		equilibrium is reached, that ad	ne voltmeter		
	13. Descri	be the function of the sa	lt bridge.		
CO	nclusions:	For the Cu/Cu ⁺² (1M), /	/ Zn ⁺² (1M), /	Zn cell: as the cell cont	inues to operate, the
ma	ss of the	electrode increase	s, and the mas	ss of the ele	ectrode decreases. The
con	centration of	ion increa	ses, and the co	oncentration of	ion decreases.